 Discovery Labs	STANDARD OPERATING PROCEDURE			
	SOP No.:	SOP-QC-013-04	Effective Date:	01.01.2017
	Supersedes:	SOP-QC-013-03	Next Review Date:	31.12.2019
	Department:	Quality Control	Page:	1 of 8
TITLE: OPERATION AND CALIBRATION TDS/CONDUCTIVITY ANALYSER				

1.0 PURPOSE:

To describe the Operation and Calibration of TDS/CONDUCTIVITY Analyser

2.0 SCOPE:

This procedure is applicable to the TDS/CONDUCTIVITY Analyzer of Laboratory.

Make : Polmon

Model : CM-540

ID No. : DIPL/QC/INS/TDS/001

3.0 RESPONSIBILITY:

3.1 Analyst-QC is responsible to follow this SOP.

3.2 Head-QC/Designee is responsible for ensuring implementation of this SOP.

3.3 Head-QA/Designee is responsible for monitoring overall compliance of this SOP.

4.0 DEFINITIONS:

Nil

5.0 PROCEDURE:

5.1 OPERATION:

5.1.1 Clean the instrument keep the area neat and clean.

5.1.2 Switch on the main power supply.

5.1.3 Press **PRG/CAL** key to enter into program mode PASSWORD menu is displayed.

PASSWORD XXX


5.1.4 Using ↑ or↓ keys to set PASS WORD (100) and use OK

Key

PRG 1. Calib

 display shows.

	Prepared by	Reviewed by	Checked by
Sign & Date			
Name	A.Navya	S.Prasad	Ch.Mahendar Reddy
Department	Quality Control	Quality Control	Quality Assurance

 Discovery Labs	STANDARD OPERATING PROCEDURE			
	SOP No.:	SOP-QC-013-04	Effective Date:	01.01.2017
	Supersedes:	SOP-QC-013-03	Next Review Date:	31.12.2019
	Department:	Quality Control	Page:	2 of 8
TITLE: OPERATION AND CALIBRATION TDS/CONDUCTIVITY ANALYSER				


- 5.1.5 To determine the **TDS of the** sample select
- PRG**
2. Option
- 5.1.6 Using ↑ or ↓ keys and use OK Key display shows. **OPTION**
- PRG**
1. READ
- 5.1.7 To select the **MEASURE MODE** press **OK** Key the display shows.
- Measure
TDS

 or

 Measure
Cond
- 5.1.8 Use ↑ Key between **TDS** or **CONDUCTIVITY**
- Measure
TDS
- 5.1.9 Select then press **OK** key
- TDS Fact
0.00
- 5.1.10 The display shows then set TDS factor 0.7 by using ↑ or ↓ keys.
- Option
2:T.Comp
- 5.1.11 Press **OK** key the display shows.
- Temp Com
ON


 Temp Com
OFF
- 5.1.12 Press **OK** Key the display using ↑ key. shows. By

	Prepared by	Reviewed by	Checked by
Sign & Date			
Name	A.Navya	S.Prasad	Ch.Mahendar Reddy
Department	Quality Control	Quality Control	Quality Assurance

 Discovery Labs	STANDARD OPERATING PROCEDURE			
	SOP No.:	SOP-QC-013-04	Effective Date:	01.01.2017
	Supersedes:	SOP-QC-013-03	Next Review Date:	31.12.2019
	Department:	Quality Control	Page:	3 of 8
TITLE: OPERATION AND CALIBRATION TDS/CONDUCTIVITY ANALYSER				

- 5.1.13 Select Temp Com
ON
- 5.1.14 Press OK key the display Temp
Auto Temp
Man shows .By using ↑ key.
- 5.1.15 Select Temp
Auto
- 5.1.16 Press OK Key the display shows then set temperature coefficient to
- T. Coeff
1% 1% by Using ↑ or ↓ or keys
- 5.1.17 Press OK key the display shows Option
3:Cell
- 5.1.18 Then press OK key the display shows Cell
1.000 adjust the cell constant given on the electrode (**Cell constant 1.02**) using ↓ or↑ keys. Then press **OK** key.
- 5.1.19 Press **PGR/CAL** key. Now instrument is ready for measuring **TDS**.
- 5.1.20 Dip the electrode/electrode assembly into **the** sample .Wait for completion analysis.
- 5.1.21 After completion the reading will display.
- 5.1.22 To determine the CONDUCTIVITY of the sample select Using ↑ or↓ key and press **OK** display shows **OPTION** menu.
- 5.1.23 To select the **MEASURE MODE** press OK key the display shows.

	Prepared by	Reviewed by	Checked by
Sign & Date			
Name	A.Navya	S.Prasad	Ch.Mahendar Reddy
Department	Quality Control	Quality Control	Quality Assurance

 Discovery Labs	STANDARD OPERATING PROCEDURE			
	SOP No.:	SOP-QC-013-04	Effective Date:	01.01.2017
	Supersedes:	SOP-QC-013-03	Next Review Date:	31.12.2019
	Department:	Quality Control	Page:	4 of 8
TITLE: OPERATION AND CALIBRATION TDS/CONDUCTIVITY ANALYSER				

**Measure
TDS**


or

**Measure
Cond**

Use ↓ or ↑ key between **TDS** or **CONDUCTIVITY**.

- 5.1.24 Select **Measure Cond** then press OK key. The display shows **Option 2:T.Comp**
- 5.1.25 Press **OK** Key the display shows. **Temp Com ON** **Temp Com OFF**
- By using ↑ key
- Select **Temp Com ON**
- 5.1.26 Press OK key the display **Temp Auto** **Temp Man** shows .By using ↑ key
- Select **Temp Auto**
- 5.1.27 Press OK Key the display **T. Coeff 1%** shows then set temperature coefficient to 1% by Using ↑ or ↓ or keys
- 5.1.28 Press OK key the display **Option 3:Cell** shows

	Prepared by	Reviewed by	Checked by
Sign & Date			
Name	A.Navya	S.Prasad	Ch.Mahendar Reddy
Department	Quality Control	Quality Control	Quality Assurance

 Discovery Labs	STANDARD OPERATING PROCEDURE			
	SOP No.:	SOP-QC-013-04	Effective Date:	01.01.2017
	Supersedes:	SOP-QC-013-03	Next Review Date:	31.12.2019
	Department:	Quality Control	Page:	5 of 8
TITLE: OPERATION AND CALIBRATION TDS/CONDUCTIVITY ANALYSER				

- 5.1.29 Then press OK key the display shows

Cell 1.000

 adjust the cell constant given on the electrode (**Cell constant 1.02**) using ↓ or ↑ keys. Then press **OK** key.
- 5.1.30 Press **PGR/CAL** key. Now instrument is ready for measuring **CONDUCTIVITY**.
- 5.1.31 Dip the electrode in to the water.

- 5.1.32 Select the AUTO RANGE by using MODE key the display

Range Auto

 shows then press OK Key
- 5.1.33 Remove the electrode from the beaker wash with distilled water and wipe with tissue paper to remove the water droplets on the surface of the electrode.
- 5.1.34 Immerse the electrode and solution whose **CONDUCTIVITY** to be measured
- 5.1.35 Record the **CONDUCTIVITY**. Value in the equipment log book.

5.2 CALIBRATION:

5.2.1 PREPARATION OF STANDARD KCL SOLUTION:

Preparation of 0.1 N KCl Solution:

Weigh accurately about 0.7492 gr of KCl into 100 ml volumetric Flask dissolve in distilled water up to the mark with the same solvent.


Preparation of 0.001 N KCl Solutions:

Pipette out exactly 1ml of 0.01N KCl solution into a 100 ml Volumetric flask dissolve in distilled water make up to the mark with the same solvent.

- 5.2.2 Press **PRG/CAL** key to enter into program mode PASSWORD menu is displayed.

PASSWORD XXX

	Prepared by	Reviewed by	Checked by
Sign & Date			
Name	A.Navya	S.Prasad	Ch.Mahendar Reddy
Department	Quality Control	Quality Control	Quality Assurance

 Discovery Labs	STANDARD OPERATING PROCEDURE			
	SOP No.:	SOP-QC-013-04	Effective Date:	01.01.2017
	Supersedes:	SOP-QC-013-03	Next Review Date:	31.12.2019
	Department:	Quality Control	Page:	6 of 8
TITLE: OPERATION AND CALIBRATION TDS/CONDUCTIVITY ANALYSER				

5.2.3 Using ↑ or ↓ keys to set PASS WORD (100) and use **OK** Key display

PRG
1. Calib

shows.

Sensor
Cal

5.2.4 To calibrate the Instrument press **OK** key the display shows

5.2.5 Press OK key the display shows scrolling message

Dip Cell in STD Solu.

5.2.6 At this stage immerse the electrode in the known conductivity solution and press

OK key the display

Wait

shows

5.2.7 After reading the stabilized **TDS/CONDUCTIVITY** segment display shows the

measured value as shows below

xxx

XX.X0C

A


Do Calib

5.2.8 The LCD display shows

5.2.9 The LED lamp(ppm/μs or ppt/ms)glows as per the sample

TDS/CONDUCTIVITY

	Prepared by	Reviewed by	Checked by
Sign & Date			
Name	A.Navya	S.Prasad	Ch.Mahendar Reddy
Department	Quality Control	Quality Control	Quality Assurance

 Discovery Labs	STANDARD OPERATING PROCEDURE			
	SOP No.:	SOP-QC-013-04	Effective Date:	01.01.2017
	Supersedes:	SOP-QC-013-03	Next Review Date:	31.12.2019
	Department:	Quality Control	Page:	7 of 8
TITLE: OPERATION AND CALIBRATION TDS/CONDUCTIVITY ANALYSER				

5.2.10 The measured value shows in segment play can adjusted to the required **STD** value as mentioned in the Table 1 by using ↑ and ↓ keys. This will automatically update cell constant.

5.2.11 Press **OK** keys to save the calibrated value. The display shows Save Y or Save N by using ↓ or ↑ keys to toggle between **Y or N**

5.2.12 Then press **OK** keys to select the option

5.2.13 Acceptance Criteria

Grams KCL/Kg of solution	K $1N\Omega^{-1}Cm^{-1}$	Acceptance Criteria
0.001N	146.9 μ	$\pm 0.5 \mu$

5.2.14 **Calibration Schedule:** Calibration shall be done every day.


6.0 FORMATS / ANNEXURE(S):

S.No	Details	Format No. (Current version)
1.	Instrument Usage log book	
2.	TDS Analyser Calibration record	

7.0 CHANGE HISTORY:

Revision No.	Effective Date	Details of Revision	Ref CCF No.
00	01.06.2007	New SOP Calibration of Equipments is introduced across all the API manufacturing facilities of Discovery.	--
01	01.08.2009	Formats changed more clear and clarity	--
02		Formats are the part of SOP. So prepared Separately.	--
03	01.01.2017		

	Prepared by	Reviewed by	Checked by
Sign & Date			
Name	A.Navya	S.Prasad	Ch.Mahendar Reddy
Department	Quality Control	Quality Control	Quality Assurance

 Discovery Labs	STANDARD OPERATING PROCEDURE			
	SOP No.:	SOP-QC-013-04	Effective Date:	01.01.2017
	Supersedes:	SOP-QC-013-03	Next Review Date:	31.12.2019
	Department:	Quality Control	Page:	8 of 8
TITLE: OPERATION AND CALIBRATION TDS/CONDUCTIVITY ANALYSER				

	Prepared by	Reviewed by	Checked by
Sign & Date			
Name	A.Navya	S.Prasad	Ch.Mahendar Reddy
Department	Quality Control	Quality Control	Quality Assurance

QA001-FM139-00